

Disposable Gowns

Certified to AAMI Class III



- 100% manufactured in the USA
- Durable, repellent, non-woven, disposable fabric
- Heavy duty nylon tie straps
- Full coverage design
- Elastic cuffs
- Standard size fits most, Oversized fits XL and above
- Fast shipment



Physical Property Summary Repellent SMS Non Woven

Property	Value Imperial	Metric	Test Method
Weight	1.26 oz/ yd ²	43.4g/m ²	WSP 130.1
Tensile	16 lbs	72N	WSP 110.1
Tear	5.0 lbs	22.5N	WSP 100.2
Penetration by Water (Hydrostatic Test)		73cm	WSP 80.6
Penetration by Water (Impact Test)		0.06g	WSP 80.3
Alcohol Repellancy	8.0 rating		WSP 80.8
Flammability	Class I		16CFR 1610.4

Model

Description

GOWN-D	Disposable Gown Standard Size Length: 44" Width: 26" (Shoulder to shoulder) Sleeve Length: 22"
GOWN-DXL	Disposable Gown Oversized Length: 44.5" Width: 27" (Shoulder to shoulder) Sleeve Length: 23.25"

CE 0338

