

Washable and Reusable Gowns

Leading Hospitals Tested & Approved!



- 100% manufactured in the USA
- Durable microfiber fabric holds up through repeated commercial washing cycles
- Can be autoclaved or ETO gas sterilized repeatedly
- Heavy duty nylon tie straps and finished edges
- Full coverage design
- Elastic cuffs
- Standard size fits most, Oversized fits XL and above
- Fast shipment



| Model | Description |
|-----------------|-----------------------------------|
| GOWN-M | Microfiber Gown Standard Size |
| | Length: 44" |
| | Width: 26" (Shoulder to shoulder) |
| | Sleeve Length: 22" |
| GOWN-MXL | Microfiber Gown Oversized |
| | Length: 44.5" |
| | Width: 27" (Shoulder to shoulder) |
| | Sleeve Length: 23.25" |