

# Revolution Maxi-Flex Loodhandschoenen



De perfecte manier om uw handen en onderarmen te beschermen door de meest flexibele, comfortabele stralingsbeschermende handschoenen op de markt.

Comfortabele pasvorm met maximale flexibiliteit. Zacht uitneembaar, wasbaar en vervangbare binnenvoeringen verminderen geuren en laten frequent wassen toe. Uitgebreide duurzame vinyl buitenhoezen zijn afneembaar voor inspectie van de loodkern en kunnen vervangen worden wanneer deze versleten zijn. De loodvinylkernen zijn nauwkeurig gevormd voor consistent lood equivalent en een lange levensduur.

Revolution Maxi-Flex loodhandschoenen komen exclusief in de Revolution-stof, het zachtste en meest comfortabele materiaal dat beschikbaar is voor stralingsbescherming.

- .25mm, .35mm, .50mm LE
- Sterke, duurzame vinyl buitenhoezen
- Zacht, flexibel en comfortabel
- Verwijderbare, wasbare en vervangbare Cool Wear vochtafvoerende binnenvoering
- Vervangbare buitenhoezen



Getoond:  
683300-705 Navy



68300-706 Teal

Code	Omschrijving
683300-701	Revolution Maxi-Flex 5 Finger Lead Gloves Red
683300-702	Revolution Maxi-Flex 5 Finger Lead Gloves Fuchsia
683300-704	Revolution Maxi-Flex 5 Finger Lead Gloves Pink
683300-705	Revolution Maxi-Flex 5 Finger Lead Gloves Navy
683300-706	Revolution Maxi-Flex 5 Finger Lead Gloves Teal
683300-700	Revolution Maxi-Flex 5 Finger Lead Gloves Black



MML | Medical

Your strategic partner in medical devices